

The Fit Test Champion

N95 DISPOSABLE RESPIRATOR



TS-H210 N95 Disposable Respirator

PURCHASE ON LINE



SCAN THIS QR CODE

Material Composition:

Filter	Polypropylene
Nose Foam	Medical Grade PVC
Nose Clip	Dual-Core Polypropylene Coated Wire
Coverweb	Polypropylene
Head Straps	Polyisoprene
Particle Filtration	98.27%
Shelf/Storage Life	3 Years from the date on the box (YYYY-MM-DD format)
Storage Temperature Limits	Between -4°F (-20°C) and +86°F (+30°C)
Storage Relativity Humidity Limits	Not Exceeding 80% RH
Country of Origin	Made in the USA
Manufactured By	Lighthouse Worldwide Solutions

Time Use Limitation:

Respirator may be used until damaged, breathing becomes difficult or contaminated with blood or bodily fluids. Discard after every use when used for surgical procedures. Follow national, state, local, and facility infection control guidance and policies.

Storage Instructions:

Store respirators in original packaging and away from contaminated areas, dust, sunlight, extreme temperatures, excessive moisture, and damaging chemicals.

Distributed By:

TISCH
Scientific



201 S. Miami Ave.
Cleveland, Ohio 44115
ScientificFilters.com, support@tisch-env.com, 513-467-0413



MADE IN USA

Lighthouse Worldwide Solutions reserves the right to change specifications without notice.